




# Microscope Camera Solutions

---

## For Evident(Olympus)Industrial Microscope

**INNOVATED FOR TOP BRAND MICROSCOPES**  
*Create a Stunning Microscope Imaging System for You*



<https://www.ostec.com.cn/>

## ■ For Trinocular Microscopes



**Applicable Models:**  
SZ61TR, SZX7, SZX10,  
SZX16, BX53M, BX53-P,  
BX51

### Option A

Measuring Camera



### Option B

WiFi Camera



### Option C

Smart Camera



### Option D

Smart Display Camera



## Function Comparison

● Standard ○ Optional – N/A

	Option A	Option B	Option C	Option D
Built-in Android OS	–	–	●	●
Pre-installed Office suits	–	–	●	●
15.6"high color gamut monitor	–	–	–	●
<b>Image output methods</b>				
5G WiFi	–	●	●	●
USB	●	●	–	–
HDMI or DP	● HDMI	● HDMI	● DP	● DP
Network	●	●	–	–

### Option A, B

- HDMI output comes with embedded software: KoPa View;
- Max 60 fps for 3840x2160 (4Kx2K) real-time preview, support snapshot and record video;
- Exclusive image modes for biological, industrial, fluorescence microscope to have accurate color reproduction;
- Option A: USB+HDMI or Network+HDMI output;  
Option B: USB+HDMI or WiFi+HDMI+Network output;
- Support depth of field fusion, stitching function (Windows software -KoPa Capture Pro);
- Support saving photos and videos to USB disk by using mouse and keyboard.

### Option C, D

- Comes with a deeply customised Android operating system, with mobile version of office suit that including Word, Excel and PPT;
- Max 30 fps for 3840x2160 (4Kx2K) real-time preview, support snapshot and record video;
- Exclusive image modes for biological, industrial, fluorescence microscope to have accurate color reproduction;
- With depth of field fusion, manual focus, stitching, measurement and other functions;
- Option C: WiF+DP output;  
Option B: local on-screen display or WiFi+DP output;
- Support saving photos and videos to USB-disk by using mouse and keyboard.

## ■ For Inverted Microscopes



Applicable Model:  
GX53

### Option A

Measuring Camera



### Option B

WiFi Camera



### Option C

Smart Camera



## Function Comparison

● Standard ○ Optional – N/A

	Option A	Option B	Option C
Built-in Android OS	–	–	●
Pre-installed Office suits	–	–	●
15.6"high color gamut monitor	–	–	–
Image output methods			
5G WiFi	–	●	●
USB	●	●	–
HDMI or DP	● HDMI	● HDMI	● DP
Network	●	●	–

## ■ For Binocular Stereo Microscopes



Applicable Models:  
SZX7, SZX10, SZX16.

Embedded beam splitter  
With 0.43X tube lens



Can be matched with the following models

### Option A

Measuring Camera  
C-mount



### Option B

WiFi Camera  
C-mount



### Option C

Smart Camera  
C-mount



## Function Comparison


● Standard ○ Optional – N/A

	Option A	Option B	Option C
Built-in Android OS	–	–	●
Pre-installed Office suits	–	–	●
15.6"high color gamut monitor	–	–	–
Image output methods			
5G WiFi	–	●	●
USB	●	●	–
HDMI or DP	● HDMI	● HDMI	● DP
Network	●	●	–


















- » HDMI output comes with measurement tools (built-in OSD menu -KoPa View);  
Support saving photos and videos to USB-disk by using mouse and keyboard.












KoPa View


















Main Functions of KoPa View

	Settings		Text		Mirror
	Color settings		Select object		Delete
	Annotation tools		Cancel		Freeze image
	Calibration		Redo		Snapshot
	Measuring tools		Flip		Record video

Camera property control

	Target brightness		White Balance		Acuity
	Automatic exposure		Contrast		Gamma
	Exposure time		Saturation		Restore default
	Gain		Chroma		

## 15 measuring tools

							
Linear distance measurement	Circular measurement	Rectangular measurement	Three-point angle measurement	Parallel line spacing measurement	Polyline measurement	Polygonal measurement	Arc measurement
							
Three-point vertical line measurement	Four-point angle measurement	Ellipse measurement	Concentric radius circle center distance drawing circle measurement	Concentric circle measurement	Ring measurement	Crosses	

## Annotation tools



Pencil



Straight line



Arrow



Rectangle




Ellipse


## » Dedicated image mode

Depending on the application of the microscope, the corresponding exclusive image mode can reproduce the image effect more accurately:

1. Biological microscope
2. Industrial microscope
3. Fluorescence microscope




### Industrial microscope



Carbon Steel (objective: 50X)



Grey Cast iron (objective: 50X)




Hypo eutectoid steel (objective: 50X)

## » Multiple image output modes

Option A: HDMI, USB, USB+HDMI, Network+HDMI Output;

Option B: HDMI, USB, WiFi, USB+HDMI, WiFi+HDMI+Network Output (USB and WiFi cannot output at the same time).





### HDMI output

Connect the camera to a monitor or large TV via HDMI cable with OSD "KoPa View".



### USB output

Connect the computer via USB /network cable by using software" KoPa Capture Pro".



### Network output

Connect the computer via network cable by using software" KoPa Capture Pro".



### 5G WiFi output(only option B)


Compatible with various devices and operating systems, including Windows, iOS, and Android. Mobile devices can access the system by scanning a QR code.

Connects to a computer via a network cable, suitable for wired long distance image transmission.



## » Depth of Field Fusion

This is achieved through Windows software: KoPa Capture Pro. Breaking through the limitation of insufficient depth of field under high-magnification objective lens, you can obtain greater depth of field by adjusting the focal length, thus obtaining sharper images than real-time single-frame images, and supporting two saves in one shot.



## Specifications

	Option A	Option B
Models	MC2000	CF48
C-mount category	AJ-C-08	BJ-C-08
With 0.43X tube lens category	AJ-A-08	BJ-A-08
Physical resolution	8.3MP	
Image sensor	SONY IMX678 CMOS	
Sensor size	1/1.8"	
Pixel size	2μm×2μm	
A/D conversion bit depth	12bit	
Exposure time	10us~10s	
Exposure mode	Rolling shutter	
ISO sensitivity	Equivalent to 100-12800	
Spectral response	400-650nm	
Exposure capability	Real-time automatic and manual adjustment	
White balance	Real-time automatic and manual RB adjustment	
Power supply	DC 5V 2A	
Video recordings	3840×2160@60fps 1920×1080@60fps	
HDMI output	3840×2160@60fps 1920×1080@60fps	
USB output	3840×2160@60fps 1920×1080@60fps	
Network output	3840×2160@60fps 1920×1080@60fps	
Software and App	Windows Software: KoPa Capture Pro, OSD: KoPa View, APP: KoPa WiFi Lab	

## Accessories

HDMI cable



USB mouse








## Dimensions(Unit:mm)

### Option A

Net weight  $\approx 0.8\text{kg}$

Camera with C-mount

Camera with 0.43X tube lens






### Option B

Net weight  $\approx 1\text{kg}$


Camera with C-mount

Camera with 0.43X tube lens




- Comes with a deeply customised Android operating system, with mobile version of office suit that including Word, Excel and PPT; support saving photos and videos to USB-disk by using mouse and keyboard.


Built-in software: KoPa WiFi Lab AO




## 18 measuring tools



- Built-in 32GB hard drive for direct saving of images and videos. Saved images can be copied via WiFi or USB-drive.



- Supports the installation of third-party office software, which can directly generate reports containing pictures and measurement results.




## Depth of field fusion

Our innovative camera leverages unique optoelectronic mechanisms and advanced image algorithms to efficiently capture images at various focal planes. This technology enables the timely synthesis of images with optimal focus, ensuring clarity and precision in every shot.


- Suitable for almost all types of microscopes on the market.
- Depth of field range, with a specific microscope, the depth of field fusion can reach max 64mm (object height).
- Fast fusion speeds and precise thresholds greatly enhance the user experience.
- During the depth fusion scanning process, the motor operates quietly and without any vibrations. The durability of the optoelectronic mechanism exceeds 3 million cycles, with each complete scan counting as one cycle.

Preview Status      Depth of field fusion results





Depth of field fusion travelling can be selected from auto- fusion, 1/3 travelling fusion, 2/3 travelling fusion, full fusion and manual fusion.

- Auto-fusion is the camera automatically scans, determines the highest and lowest points on the object side, and fuses every clear point of the fallout.
- 1/3 travelling fusion is a forced scan fusion of one third of the camera's physical maximum scan stroke.
- 2/3 traveling fusion is a forced scanning fusion of two thirds of the camera's physical maximum scanning stroke.
- Full fusion is a forced scanning fusion for the physical maximum scanning stroke of the camera.
- Manual fusion is a scanning fusion of the region of interest by controlling the scanning stroke with the mouse wheel.



## » Manual focus


Manually roll the mouse wheel and click on the area you want to observe to focus on it.




## » Dedicated image mode

Depending on the application of the microscope, the corresponding exclusive image mode can reproduce the image effect more accurately:


1. Biological microscope
2. Industrial microscope
3. Fluorescence microscope




### Industrial microscope



Carbon Steel (objective: 50X)



Grey Cast iron (objective: 50X)





Hypo eutectoid steel (objective: 50X)

## » Multiple image output modes

Option C: DP, WiFi, WiFi+DP output;

Option D: local on-screen display, DP, WiFi, WiFi+DP output.



### 5G WiFi output

Compatible with various devices and operating systems, including Windows, iOS, and Android. Mobile devices can access the system by scanning a QR code. connects to PC via WiFi.




### DP output

With DP output for display on monitors, TVs and projectors.



### Connecting to the Internet

The camera supports wireless Internet connection (only supports 5G WiFi router signal), by entering the password of any 5G WiFi network, you can directly open the browser to access the Internet. This will allow you to take advantage of online features and functionalities directly from your camera.





## Specifications

### 15.6" high color gamut display

Number of pixels	1920( horizontal) x 1080 (vertical)
Pixels arrangement	RGB vertical stripe
Colour gamut	100% (sRGB)
Display number of colors	16.7M(8Bit)
Surface treatment	Anti-glare
Surface hardness	3H
Viewing angel range	170 horizontal, 170 vertical
Contrast	800
Brightness	500cd/m <sup>2</sup> (average of 5 points)

	Option C	Option D
Models	TE2000	JX200
C-mount category	CJ-C-08	DJ-C-08
With 0.43X tube lens category	CJ-A-08	DJ-A-08
Physical resolution	8.3MP	
Image sensor	SONY IMX678 CMOS	
Sensor size	1/1.8"	
Pixel size	2μm×2μm	
A/D conversion bit depth	12bit	
Exposure time	10us~10s	
Exposure mode	Rolling shutter	
ISO sensitivity	Equivalent to 100-12800	
Spectral response	400-650nm	
Exposure capability	Real-time automatic and manual adjustment	
White balance	Real-time automatic and manual RB adjustment	
Power supply	DC 12V 3A	
Video recordings	3840×2160@30fps 1920×1080@30fps	
HDMI output	3840×2160@30fps 1920×1080@30fps	
USB output	3840×2160@30fps 1920×1080@30fps	
Network output	3840×2160@30fps 1920×1080@30fps	
Software and App	Windows Software:KoPa Capture Pro,Embedded software:KoPa WiFi Lab AO,App:KoPa WiFi Lab	

## Accessories

DP cable  
(option C)



Power adapter and power cord  
(Optional Chinese, American, European,  
Australian, Korean, British standard etc.)




USB mouse and  
keyboard




## Dimensions(Unit:mm)

### Option C


Net weight  $\approx 1.3\text{kg}$



Camera with 0.43X tube lens




Camera with C-mount




### Option D

Net weight  $\approx 2\text{kg}$

Camera with 0.43X tube lens



Camera with C-mount



1. Comply with FCC certification of The US Federal Communication Commission.
2. Comply with European (standard) safety CE certification.
3. Comply with the MIC certification issued by the Ministry of Internal Affairs and Communications of Japan (Electric Wave Method and Electro-Optical Communication Business Law).
4. Comply with JATE certification of Japanese telecommunications law directive.
5. Comply with the "Directive on the Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment" (RoHS) Directives in accordance with EU legislation.

Evaluation object	Certification	Certificate File Name & Report	Certificate number & corresponding report number
WF01A(5G WiFi 11ac)module Certification	US FCC Report	SZEM180100024801-5G wifi RPT-WF01A FCC Report	SZEM180100024801
		SZEM180100024802-RT-WF01A FCC Report	SZEM180100024802
		Appendix A-Photographs of EUT Constructional Details for SZEM1801000248CR-FCC	SZEM1801000248CR
	US FCC ID Certification	2AFO3WF01A_NII-WF01A FCC ID	2AFO3WF01A
	EU CE report	SZEM180100024901 EN301489 RPT-WF01A CE Report	SZEM180100024901
		SZEM180100024902 WIFI5G RPT-WF01A CE Report	SZEM180100024902
	Japanese MIC Certification	CSRT180084-WF01A Japanese MIC Certification	CSRT180084
	Japanese JATE Certification	CSTT180018-WF01A Japanese JATE Certification	CSTT180018

## Patented

Patent category	Patent name	Patent number
Design patent	Electronic eyepiece	ZL 2015 3 0193227.8
	Wireless electronic eyepiece	ZL 2015 3 0193223.X
	Electronic eyepiece with spectroscopic system	ZL 2019 3 0331144.9
	Microscope (with splitting prism camera)	ZL 2019 3 0717439.X
	Microscope with camera	ZL 2019 3 0717442.1
Utility model patents	WiFi microscope eyepiece	ZL 2015 2 0296469.4
	Electronic eyepiece	ZL 2015 2 0426409.X
	Wireless electronic eyepiece	ZL 2015 2 0426313.3
	Microscope with displayer	ZL 2019 2 0928962.1
	Electronic eyepiece with splitting prism system	ZL 2019 2 1022863.3

## Software copyright

Category	Name of software	Platform	License number
Computer software copyright registration certificate	KoPa Capture Pro	Windows	2021SR1287730
	KoPa WiFi Lab AO	Android	2021SR1304520
	KoPa WiFi Lab	Android	2019SR0117768
		iOS	2019SR0028558
	KoPa View	Linux	2024SR1617066

**KoPa®** GuangZhou Ostec Electronic Technology Co.,Limited

Manufacturer: No.8 West Lane, Jiangcheng Road, Bangjiang East Village,Dalong street, Panyu District, Guangzhou, China.



High-Tech Enterprise certificate number:  
GR202344009665



ISO9001 Verification No:00223Q26818R3S

The content of this leaflet has been reviewed by our company at the time of its release. Due to technological development, the actual product is subject to change without notice.

The names of other companies, product names, and trademarks **OLYMPUS** **Nikon** **Leica** **ZEISS** **Apple** **Android** **HarmonyOS** **W** **Q** **du** recorded on this leaflet are owned by their companies